

GMFCS  
Levels III-IV

kg  
max. 100



46-124 cm

29-113 cm



### I am your versatile posterior walking

**aid.** I am always there for you - even from your first year of life. Walking is my biggest hobby. It doesn't matter where you want to go; I will follow you. If you support yourself on my handles, you will quickly notice how I take on your body weight. In addition, my special abducted frame shape supports you so that you don't swing outwards when walking. My pelvic support pad ensures that you walk upright and straight. Cool, aren't I?!





**Combined pelvic and buttock pad**  
can become a backrest in one  
single easy step.

**Forearm rests**  
with handles.

**Padded seat**  
can easily be folded back.

**Aluminium frame**  
can easily be folded up.

**Direction locks**  
facilitate  
walking straight.

**Backwards stop**  
prevents inadvertent  
rolling back.

**EVA guide wheels**  
make the walking aid manoeuvrable.



**Frame colours**  
Size 0 mini



**Frame colours**  
Size 1



**Frame colours**  
Size 2



**Frame colours**  
Size 3 and Size 4





Data	Size 0 mini*	Size 1	Size 2	Size 3	Size 4
Handle height	29 - 53 cm	34 - 58 cm	51 - 75 cm	68 - 94 cm	87 - 113 cm
Height of forearm rests	-	46 - 70 cm	63 - 87 cm	79 - 103 cm	98 - 124 cm
Handle spacing	22 - 36 cm	22 - 36 cm	30 - 44 cm	38 - 52 cm	43 - 57 cm
Seat height	-	33 cm	45 cm	60 cm	70 cm
Total width	52 cm	51 cm	61 cm	68 cm	75 cm
Total length	58 cm	58 cm	70 cm	77 cm	87 cm
Frame height (folded)	37 cm	25 cm	32 cm	34 cm	37 cm
Transport dimensions (L x W x H)	58 x 52 x 37 cm	58 x 51 x 25 cm	70 x 61 x 32 cm	77 x 68 x 34 cm	87 x 75 x 37 cm
max. load	20 kg	30 kg	50 kg	75 kg / 70 kg	100 kg
Weight	3.9 kg	7.1 kg	8.3 kg	8.5 kg / 9.8 kg	10.2 kg
HMV No.	10.46.02.3029				

\* Frame cannot be folded up



You can find more information on accessories in the price list.

